printo[®] dent Generative Resin GR-23 mouthguard

Instructions for use
Gebrauchsanweisung
Mode d'emploi
Instrucciones de uso
Istruzioni per l'uso

pro^{3d}ure
medical
Materials are our DNA!

1. Intended use / Indication

Light-curable polymerizable resin intended to be used in conjunction with extra-oral curing light equipment for the fabrication, by additive manufacturing, of dental appliances. Indicated for the fabrication of mouthguards.

2. Contraindication

printo[®] GR-23 mouthguard is contraindicated ...

- ... if a patient is known to be allergic to any of the ingredients.
- ... for every application that is not part of the indication (see above).

3. Patient Target Group

Sportspeople being treated in the context of a dental procedure.

4. Intended Users

Dentist, dental technicians

5. Requirements

pro^{3d}ure CL-1 is mandatory for post cleaning process!

Software - information available from:

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Hardware (3D Printing) - information available from:

ASIGA · 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia

Hardware (Post Curing) - information available from:

NK Optik GmbH (otoflash G171) · Isarstr. 2, D-82065 Baierbrunn, Germany

Additional specifications on www.pro3dure.com.

6. Material

printo[®] GR-23 mouthguard consists of functional (meth)acrylic resins, initiators, dyes and stabilisers.

7. Geometric presetting

Minimum wall thickness: 2 mm

8. Material parameters

Validated parameters for the print of printo[®] GR-23 mouthguard are provided as an INI-file in the ASIGA composer software on the ASIGA MAX UV printer. Please make sure to select the right file for your product.

9. Manufacturing process (fig. 1-13)

- Prepare data (CAD & build preparation).
- Choose process parameters (build style, etc.).
- Transfer prepared data to 3D printing device.
- Prepare 3D printing – shake bottle.
- Fill resin tank of 3D printing device.
- Build the parts.
- Remove most of the excess resin from the green parts using a cloth.
- Sonicate the green prints in CL-1 for 5 minutes.
- Brush with CL-1 any areas where the sonication did not remove the resin.
- Then sonicate in IPA for 5 minutes.
- Rinse with IPA if necessary.
Note: IPA sonication is only required to remove CL-1.
- Post-cure in Otofash G 171: 2 x 4000 flashes in inert atmosphere.
- Final cleaning: 30 minutes water sonication.

10. Finishing processes

Polishing


11. Notice

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred. **Caution:** The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not handle until all safety precautions have been read and understood. Avoid breathing dust/fume/gas/mist/vapours/spray. Do not eat, drink or smoke when using this product. Use only outdoors or in a well-ventilated area. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. IF ON SKIN: Wash with plenty of soap and water. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Rinse mouth. If skin irritation occurs: Get medical advice/attention. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Collect spillage. Store in a well-ventilated place. Keep container tightly closed. Store locked up. Dispose of contents/container according to official regulations.

12. Hazard statements

Harmful if swallowed. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye damage. May cause respiratory irritation. Very toxic to aquatic life with long lasting effects.

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
Manufacturer:
pro^{3d}ure medical GmbH
Am Burgberg 13 · 58642 Iserlohn, Germany
Phone: +49 (0)2374 920050-0

Distributor (US):
pro^{3d}ure medical LLC
9825 Valley View Road · Eden Prairie, MN 55344
Phone: 952-426-1928

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